





CARDIOVASCULAR
Embolectomy catheters



Value Life

Embolectomy catheters

- Short and flexible end allowing easy introduction into the artery while reducing the risks of plaque lesion
- Strong balloon with cylinder shape when inflated



Features & benefits

BALLOON IN LATEX

COLOUR CODED CATHETER

- For easy identification of the French sizes

PRINTED CATHETER

- Maximum balloon inflation capacity printed on the proximal end of the catheter for safety
- Graduation every 10cm

References

| References | Size (Fr) | Length (cm) | Colour | Max. fluid capacity (ml) | Inflated balloon diameter (mm) |
|------------|-----------|-------------|--------|--------------------------|--------------------------------|
| CEM402 | 2 | 40 | | 0.20 (Gaz) | 4 |
| CEM602 | 2 | 60 | | 0.20 (Gaz) | 4 |
| CEM802 | 2 | 80 | | 0.20 (Gaz) | 4 |
| CEM403 | 3 | 40 | | 0.20 | 5 |
| CEM603 | 3 | 60 | | 0.20 | 5 |
| CEM803 | 3 | 80 | | 0.20 | 5 |
| CEM404 | 4 | 40 | | 0.75 | 7 |
| CEM604 | 4 | 60 | | 0.75 | 7 |
| CEM804 | 4 | 80 | | 0.75 | 7 |
| CEM405 | 5 | 40 | | 1.50 | 9 |
| CEM605 | 5 | 60 | | 1.50 | 9 |
| CEM805 | 5 | 80 | | 1.50 | 9 |
| CEM406 | 6 | 40 | | 2.00 | 11 |
| CEM806 | 6 | 80 | | 2.00 | 11 |
| CEM407 | 7 | 40 | | 2.50 | 13 |
| CEM807 | 7 | 80 | | 2.50 | 13 |
| CEM808 | 8 | 80 | | 3.00 | 14 |

Please read carefully the Instructions For Use. The CE conformity assessment has been conducted by GMED (0459).

 **CARDIOVASCULAR & SURGERY**

For further information, please contact: questions@vygon.com

The specifications shown in this leaflet are for information only
and are not, under any circumstances, of a contractual nature.

Vygon – 5, rue Adeline • 95440 ECOUEN • FRANCE
Reception: +33 (0)1.39.92.63.63 – Service clients France: +33 (0)1.39.92.63.81
Export customer service: +33 (0)1.39.92.64.15
Fax.: +33 (0)1.39.92.64.44 • www.vygon.com

